



THE STERISHIELD DELIVERY SYSTEM

The closed and validated trigger spray system that ensures the integrity of the contents throughout use

INTRODUCING THE STERISHIELD DELIVERY SYSTEM

The Ideal Solution for Cleaning and **Disinfection of Critical Areas**



An effective cleaning and disinfection programme in controlled areas is critical to maintain product quality and safety in the pharmaceutical manufacturing industry.

Typical delivery solutions such as trigger spray systems and aerosols can compromise various areas of your business, including:

Employee Health and Safety

Cost of Use / Profitability

Product Quality



Regulatory Compliance

"Disinfectants should be shown to be effective for the duration of their in use shelf-life."

STERISHIELD DELIVERY SYSTEM (SDS)

Specifically designed to address cleaning and disinfection challenges, the SteriShield Delivery System's unique combination of the patented trigger system and integral bag provides a closed system and ensures contamination cannot enter the bottle during use. Unique features of this revolutionary design include:

- Closed System No air suck back ensures sterility
- Validated System Sterile in-use shelf life
- ▲ Variable Spray Pattern Ensures optimum surface coverage
- ▲ Reduced Airborne Particulates Safer compared to aerosols





AEROSOLS

TRIGGER SYSTEMS

Ê

Ŧ

Do not form a closed system and are unable to prevent air from being drawn back into the bottle, posing a risk of contaminating the contents of the bottle during use. This is applicable to both bagin-bottle and standard trigger spray systems. This puts your product quality and compliance at risk due to potential loss of sterility of the contents and increases in-use cost.

Can pose a safety risk due to the generation of higher levels of small airborne particulates and may result in an increase in overall cost given their higher disposal

costs. Aerosols also lack an adjustable spray pattern that may limit the coverage of surface areas.

WHAT MAKES THE TRIGGER SO EFFECTIVE?

A patented precision engineered trigger head to ensure that air is not drawn back into the bottle

WHAT MAKES THE BOTTLE SO **EFFECTIVE?**

A bottle with dual layers:

- ▲ A rigid outer bottle to provide structural stability
- ▲ An integral bag protects its contents and creates a postdispensing vacuum in which remaining liquid is stored and kept sterile

Reference: Eudralex, new Annex 1 (draft), 5.31

PHARMACEUTICAL



C

FEATURES

A variable rotating nozzle to deliver the desired spray pattern

A shroud that is smooth, simple and easy to disinfect

An ergonomically designed trigger to provide ease of use and consistency for every operator

Optimum dosing volume

Precision engineered pump system, lubricated with medical grade silicone, to eliminate air ingress around the plungers

A unique venting system, fitted with a 'duck billed valve', to ensure ONLY one way venting during irradiation

An inner bag that is sealed to the outer bottle all around the top of the bottle neck

An inner flexible bag that is co-extruded with the bottle when manufactured

A plastic outer bottle which provides structural rigidity

Two holes drilled into the base of the outer bottle to allow air to move freely into the space between the flexible bag and the outer bottle, enabling the inner bag to collapse around the fluid

An inner bag is sealed to the outer bottle at the centre of the base to allow it to collapse around the dip tube

Why use the Ecolab SteriShield **Delivery System** Trigger Spray?

The SteriShield Delivery System trigger spray delivers the benefits of sterility by preserving the integrity of the product, minimising risk to operators, fulfilling regulatory requirements and providing a cost-effective solution.



HEALTH & SAFETY

- ▲ The variable nozzle allows for a larger spray droplet size to reduce atmospheric concentration of products such as alcohol
- No pressurised content



COST OF USE

- ▲ Uses the entire volume of product
- Product can be used for the entire in-use shelf life



PRODUCT QUALITY

- ▲ Closed system eliminating contamination risk
- Adjustable spray pattern ensuring optimum surface coverage
- Easy to clean shroud and bottle



REGULATORY COMPLIANCE

- Validated in-use content sterility and shelf life
- Validated trigger spray



EFFICACY ASSURED

The variable nozzle allows the liquid to be dispensed as required and the correct application can be obtained for each procedure enabling thorough wetting of the surface so that all disinfection procedures are effective.

TESTED AND VALIDATED

The SDS has been subjected to a stringent and comprehensive validation testing programme. The system has undergone in-depth validation, including microbiology media broth tests to detect any potential ingress of viable contamination from the external environment.

Additionally, Ecolab has undertaken vacuum and particulate testing to prove without doubt that the trigger and bottle combination creates and maintains a closed system, protecting the sterility of the contents in-use.

OPERATOR SAFETY

The SDS offers improved safety compared to aerosol sprays. The larger droplet size dispensed from the trigger spray reduces the level of atmospheric alcohol created at the time of use, consequently reducing the risk of excessive vapour inhalation by operators, while also reducing false positives in air sampling.

PHARMACEUTICAL

OPTIMUM SURFACE COVERAGE

The fully adjustable trigger spray nozzle allows the liquid to be dispensed in any configuration, from a jet to a wide spray, allowing the user to select the most appropriate method of application.

REDUCED COST OF USE

By providing a closed system, the SDS offers reduced cost by eliminating expensive, pressurised containers with their associated waste disposal requirements and wastage of unused product which cannot be dispensed.

The system also removes the risk of contamination associated with other trigger systems.

From purchase to disposal, the SDS can offer a significant advantage by reducing the total cost of use.

EASY TO USE

The development of our user friendly trigger provides:

- ▲ The same quality product from start to finish
- An adjustable nozzle
- ▲ Variable spray pattern for specific requirements
- An easy to clean shroud
- A Negligible wastage
- Ergonomic design ensuring comfortable use for a range of hand sizes



Always look for the SDS stamp on our packaging the patented and validated trigger spray system

Validation Studies for the **SteriShield Delivery System**



Note: The successful validation test results for the SDS are not applicable to any other trigger spray system.

How do you prove the SteriShield System is a closed system?

As part of our development of the SteriShield Delivery System, the system was subjected to rigorous testing to ensure its performance as a closed system.

A combination of vacuum tests, nutrient broth tests and particle tests have been conducted with results summarised below. A full set of data is available for each of these tests.

The results provide the assurance of a fully validated dispensing system which is fitted to our spray products - ensuring long in-use shelf life and peace of mind. "The combination of Ecolab's unique trigger head and bottle are what makes the SDS a truly closed system."



NUTRIENT BROTH TEST

Bottles are filled with microorganism growth promotion media (both FTM and TSB), periodically used (sprayed) over a six month period and the contents regularly tested for continued sterility. In total 90 bottles were tested to ensure robust data. All test samples showed no growth, indicating sterility

PARTICLE TEST

An empty SDS bottle with a trigger spray is connected to a particle counter and the trigger is 'pumped' multiple times to simulate use.

The system measures 5 and 0.5 micron particles.

Various leading competitor products were subjected to the same particle test and all **failed** (as shown in the table bottom right).

- The system was self validating (as 0 was required before each test)
- The test is a demonstration of real time in-use situations

The results highlight the risk of not using the SDS



Particle Count Zero

Pump

TRIGGER

SDS Trigger Spray

Competitor 1

Competitor 2

Competitor 3

*In these two test



A trigger spray without the SDS showed

growth with 1/10th of volume usage.



> 0.5 microns particles			> 5 microns particles		
Test 1	Test 2	Test 3	Test 1	Test 2	Test 3
0	0	0	0	0	0
0	0	0	0	0	0
0	0	0	0	0	0
0	0	0	0	0	0
421	15	103	4.2	1.1	0
0	0	0	0	0	0
249	64	123	0	2.1	2.1
368	82	82	24.6	0	0
105	100	150	2.1	0	0
28	26	74*	4	4	15*
62	353	56	0	21	8
189	25	251	21	1	29

*In these two tests it was noted that the trigger and bottle combination were unable to clean up after use. This indicated that the use of the trigger had compromised it permanently. (There was no evidence of any damage to the trigger.)

7



Our programme consists of a wide array of products and services that meet your key needs for cleanroom cleaning and disinfection to ensure patient safety, product quality, regulatory compliance, employee health and safety and operational efficiency.

PLEASE SPEAK TO YOUR ECOLAB ACCOUNT MANAGER FOR FURTHER INFORMATION

WORLDWIDE HEADQUARTERS 1 Ecolab Place St. Paul, MN 55102 USA www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7 8304 Wallisellen Switzerland www.ecolablifesciences.com

